

# SMALL MAMMAL ENDOSURGERY



Dr Laila Proença

UGA -Supreme Pet Foods Resident in Zoological Medicine  
MV, DVM, Ms, PhD



The University of Georgia





1st International Conference on Avian, Herpetological and Exotic Mammal Medicine 2013  
(12th European AAV Conference, 2nd International ARAV Conference, 1st International AEMV Conference and 2nd ECZM Scientific Meeting)

April 20-26, 2013 - Wiesbaden, Germany









# INDICATIONS




- Ovarian neoplasia, uterine neoplasia, pancreas, lymph
- Gastric neoplasia, compared with
- Ovarian neoplasia, oophorectomy
- Tumor resection, gastropexy, cystopexy, urinary calculi removal, adrenalectomy, insulinoma resection, laparoscopic-assisted GI surgery



# INDICATIONS






- Ovarian neoplasia, uterine neoplasia, pancreas, lymph
- Gastric neoplasia, compared with
- Ovarian neoplasia, oophorectomy
- Tumor resection, gastropexy, cystopexy, urinary calculi removal, adrenalectomy, insulinoma resection, laparoscopic-assisted GI surgery





# INDICATIONS





neoplasias, lymph

associated with

stomach

- Gastropexy, cystopexy, urinary calculi removal, adrenalectomy, insulinoma resection, laparoscopic-assisted GI surgery

# INDICATIONS






- Ovarian cysts, uterine masses, lymph
- Gastric masses, compared with
- Ovarian cysts, oophorectomy
- T
- Gastropexy, cystopexy, urinary calculi removal, adrenalectomy, insulinoma resection, laparoscopic-assisted GI surgery






# CONTRAINDICATIONS




- Same as contraindications for laparotomy
- Lack of appropriate endoscopy equipment
- Lack of appropriate training and experience

# PATIENT POSITIONING



# PATIENT POSITIONING


- Aim to present the surgical site at the highest level
  - Raising the head (reverse Trendelenburg): Helps with cranial surgery (e.g. liver evaluation)






# PATIENT POSITIONING

- Aim to present the surgical site at the highest level
  - Raising the head (reverse Trendelenburg): Helps with cranial surgery (e.g. liver evaluation)






# PATIENT POSITIONING

- Aim to present the surgical site at the highest level
  - Raising the head (reverse Trendelenburg): Helps with cranial surgery (e.g. liver evaluation)
- Lowering the head (Trendelenburg)
  - Helps with caudal surgery (e.g. bladder, retained testis)








# PATIENT POSITIONING



- Dorsal recumbency for ventral structures
  - e.g. liver, bladder, stomach
- Dorsolateral to lateral recumbency for dorsolateral structures
  - e.g. ovary, kidney, adrenal
- Tilting tables
- Positioners
- Towels and tape




# PATIENT POSITIONING






# TELESCOPE-CANNULA ORIENTATION

- Patient preparation
- Classic base-ball infield orientation
  - Promotes triangulation
  - Avoids parallel instruments
  - Provides some depth perception
- Experienced endosurgeons can violate this principle
  - Tedious, fatiguing and frustrating!



# TROCAR-CANNULA PLACEMENT

- **Two techniques for initial insufflation**
- Veress needle
  - Fast
  - Greater danger of GI perforation
- Hasson Technique
  - Slower and requires suture
  - Safer placement






# TROCAR-CANNULA

- **Hasson Technique**

- 3 mm cannula and blunt dissector
- 1 cm mini-laparotomy on ventral midline
- Skin and SQ fat incised
- Traction sutures placed and used to reduce leakage
- Elevation of linea alba and small incision
- Use curved forceps to ensure peritoneal perforation
- Replace with cannula, forced to right cranial quadrant




can be

# TROCAR-CANNULA

- **Hasson Technique**

- 3 mm cannula and b
- 1 cm mini-laparotomy  
ventral midline
- Skin and SQ fat in
- Traction sutures p  
used to reduce lea
- Elevation of linea alba and small incision
- Use curved forceps to ensure peritoneal perforation
- Replace with cannula, forced to right cranial quadrant






can be



# TROCAR-CANNULA


- **Hasson Technique**

- 3 mm cannula and 1 cm mini-laparotomy
- 1 cm mini-laparotomy ventral midline
- Skin and SQ fat incision
- Traction sutures used to reduce leakage
- Elevation of linea alba and small incision
- Use curved forceps to ensure peritoneal perforation
- Replace with cannula, forced to right cranial quadrant



# TROCAR-CANNULA PLACEMENT

- **Hasson Technique**
  - Place purse-string or mattress suture around cannula to improve security and air-tight seal
- Midline only






# TROCAR- CANNULA PLACEMENT


- Second and third cannulae

- Craniocaudal along midline
- Useful for lateral structures (e.g. oophorectomy)




- Second and third cannulae

- Transverse through muscles
- Useful for cranial and caudal procedures (e.g. liver biopsy)



# TROCAR-CANNULA PLACEMENT

- Initial insufflation pressure set to 8-10 (12) mmHg
  - 1-2 L/min for rabbits
  - 0.1-0.5 L/min for rodents
- Creates a turgid abdominal wall that facilitates 2nd and 3rd trocar-cannulae placement under telescope guidance
  - Transillumination of skin/muscle helps avoid vessels






# INSUFFLATION

- Once cannulae are in place reduce insufflation pressure to 8 mmHg
- Flow rates can also be reduced but are best controlled automatically using electronic endoflator



# RABBIT LAPAROSCOPY - SINGLE PORT





# RABBIT LAPAROSCOPY - DOUBLE PORT



## Rabbit Laparotomy Dual instrument liver biopsy



Prof Stephen Divers & Dr Laila Proenca


# OOPHORECTOMY

- Three port technique
- Dorsolateral position
- Grasp ovary using inferior instrument (atraumatic forceps) and elevate
- Use bipolar forceps or monopolar scissors (superior instrument) to coagulate vascular supply
- Use monopolar scissors to dissect ovary free








# OOPHORECTOMY



# CLOSURE

- Evacuation of all gas
- Skin closure
  - single suture
  - tissue adhesive
- Suture removal does not result in dehiscence or evisceration.








# OUTCOMES

- **Advantages**

- Less invasive
- Less painful
- Faster return to normal behaviors




- **Disadvantages**

- Requires more equipment
- Endosurgery often requires a 2nd surgeon or assistant



# TRAINING & COMPETENCE




- **Endosurgery courses**
  - [www.vet.uga.edu/mis](http://www.vet.uga.edu/mis)
- **Literature**
  - Small Animal Endoscopy (Jan 2011)
  - Vet Clinic North America (May 2010)
- **Laparoscopy trainers**
  - Home-made
- **Give yourself options!**
  - Start with endosurgery but always retain the option to convert to traditional procedure!





# TRAINING & COMPETENCE

- Endosurgery courses
  - [www.vet.uga.edu/mis](http://www.vet.uga.edu/mis)
- Literature
  - Small Animal Endoscopy (Jan 2011)
  - Vet Clinic North America (May 2010)
- Laparoscopy trainers
  - Home-made
- Give yourself options!
  - Start with endosurgery but always retain the option to convert to traditional procedure!





# IN CONCLUSION...

- Laparoscopy offers a minimally invasive approach to the abdomen of small mammals
- Facilitates biopsy, remote injection techniques, tissue removal/resection
- With practice, laparoscopy times are shorter than laparotomy





# ACKNOWLEDGEMENT

- Supreme Pet Foods
- UGA
- Dr Stephen Divers
- Dr Jörg Mayer
- ICARE
- Karl Storz

